

# NIFK Antibody



PACO05516

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

NIFK

**Uniprot**

Q9BYG3

**Synonyms:**

MKI67IP; NIFK; NOPP34; MKI67 FHA domain-interacting nucleolar phosphoprotein; Nucleolar phosphoprotein Nopp34; Nucleolar protein interacting with the FHA domain of pKi-67; hNIFK

**Immunogen:**

Synthesized peptide derived from human NIFK around the non-phosphorylation site of T234.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A