

GPR174 Antibody



PACO05530

Product Information

| | |
|--|--|
| Size: 50ug | Protein Background: |
| Reactivity: Human | Gene ID: GPR174 |
| Source: Rabbit | Uniprot Q9BXC1 |
| Isotype: IgG | Synonyms: GPR174; FKSG79; Probable G-protein coupled receptor 174 |
| Applications: ELISA, WB, IF | Immunogen: Synthesized peptide derived from the Internal region of human GPR174. |
| Recommended dilutions: ELISA:1:40000, WB:1:500-1:2000, IF:1:200-1:1000 | Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |

Product Images

N/A

N/A