

RERG Antibody



PACO05553

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

RERG

Uniprot

Q96A58

Synonyms:

RERG; Ras-related and estrogen-regulated growth inhibitor

Immunogen:

Synthesized peptide derived from the Internal region of human RERG.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A