

TCEAL3/TCEAL5/TCEAL6 Antibody



PACO05555

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TCEAL3/TCEAL5/TCEAL6

Uniprot

Q969E4/Q5H9L2/Q6IPX3

Synonyms:

TCEAL3; MSTP072; Transcription elongation factor A protein-like 3; TCEA-like protein 3; Transcription elongation factor S-II protein-like 3; TCEAL5; Transcription elongation factor A protein-like 5; TCEA-like protein 5; Transcription elonga

Immunogen:

Synthesized peptide derived from the C-terminal region of human TCEAL3/5/6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A