

Phospho-IKBKG (S85) Antibody



PACO05559

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

IKBKG

Uniprot

Q9Y6K9

Synonyms:

IKBKG; FIP3; NEMO; NF- κ B essential modulator; NEMO; FIP-3; I κ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor κ B kinase subunit gamma; I- κ B kinase subunit gamma; IKK-gamma; IKKG; I κ B kinase subunit gamma; NF

Immunogen:

Synthesized peptide derived from human IKK γ ; around the phosphorylation site of S85.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A