

# HAT1 Antibody



PACO05561

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

HAT1

**Uniprot**

O14929

**Synonyms:**

HAT1; KAT1; Histone acetyltransferase type B catalytic subunit; Histone acetyltransferase 1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human HAT1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A