

PACO05562

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

DGKZ

**Uniprot**

Q13574

**Synonyms:**

DGKZ; DAGK6; Diacylglycerol kinase zeta; DAG kinase zeta; Diglyceride kinase zeta; DGK-zeta

**Immunogen:**

Synthesized peptide derived from the Internal region of human DGK- zeta.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A