

BAGE2 Antibody



PACO05566

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:

Gene ID:

BAGE2

Uniprot

Q86Y30

Synonyms:

BAGE2; B melanoma antigen 2; Cancer/testis antigen 2.2; CT2.2

Immunogen:

Synthesized peptide derived from the Internal region of human BAGE2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A