

# OR6A2 Antibody



PACO05583

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR6A2

**Uniprot**

O95222/Q3MJC7

**Synonyms:**

Olfactory receptor 6A2

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 6A2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A