

DDX3Y Antibody



PACO05595

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

DDX3Y

Uniprot

O15523

Synonyms:

DDX3Y; DBY; ATP-dependent RNA helicase DDX3Y; DEAD box protein 3; Y-chromosomal

Immunogen:

Synthesized peptide derived from the N-terminal region of human DDX3Y.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A