

B4GALT2 Antibody



PACO05610

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

B4GALT2

Uniprot

O60909

Synonyms:

B4GALT2; beta-1; 4-galactosyltransferase 2; beta-1,4-GalTase 2; beta4Gal-T2; b4Gal-T2; UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 2; UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 2

Immunogen:

Synthesized peptide derived from the C-terminal region of human β -1,4-Gal-T2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A