

B3GALT2 Antibody



PACO05611

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

B3GALT2

Uniprot

O43825

Synonyms:

B3GALT2; beta-1; 3-galactosyltransferase 2; beta-1,3-GalTase 2; beta3Gal-T2; beta3GalT2; UDP-galactose:2-acetamido-2-deoxy-D-glucose 3beta-galactosyltransferase 2

Immunogen:

Synthesized peptide derived from the C-terminal region of human β -1,3-Gal-T2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A