

TNFSF14 Antibody



PACO05618

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:20000, IF:1:200-1:1000

Protein Background:

Gene ID:

TNFSF14

Uniprot

O43557

Synonyms:

TNFSF14; HVEM-L; LIGHT; Tumor necrosis factor ligand superfamily member 14; Herpes virus entry mediator ligand; HVEM-L; Herpesvirus entry mediator ligand; CD antigen CD258

Immunogen:

Synthesized peptide derived from the Internal region of human LIGHT.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A