

FGF18 Antibody



PACO05626

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

FGF18

Uniprot

O76093

Synonyms:

FGF18; Fibroblast growth factor 18; FGF-18; zFGF5

Immunogen:

Synthesized peptide derived from the C-terminal region of human FGF-18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A