

CCKAR Antibody



PACO05635

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CCKAR

Uniprot

P32238

Synonyms:

CCKAR; CCKRA; Cholecystokinin receptor type A; CCK-A receptor; CCK-AR;
Cholecystokinin-1 receptor; CCK1-R

Immunogen:

Synthesized peptide derived from the Internal region of human CCK-AR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A