

CEACAM21 Antibody



PACO05667

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CEACAM21

Uniprot

Q3KPI0

Synonyms:

CEACAM21; Carcinoembryonic antigen-related cell adhesion molecule 21

Immunogen:

Synthesized peptide derived from the Internal region of human CEACAM21.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A