

Phospho-DOK2 (Y345) Antibody



PACO05671

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

DOK2

Uniprot

O60496

Synonyms:

DOK2; Docking protein 2; Downstream of tyrosine kinase 2; p56(dok-2)

Immunogen:

Synthesized peptide derived from human Dok-2 around the phosphorylation site of Y345.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A