

CLDN10 Antibody



PACO05679

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CLDN10

Uniprot

P78369

Synonyms:

CLDN10; Claudin-10; OSP-like protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human Claudin-10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A