

Phospho-CLDN6 (Y219) Antibody



PACO05681

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CLDN6

Uniprot

P56747

Synonyms:

CLDN6; Claudin-6; Skullin

Immunogen:

Synthesized peptide derived from human Claudin-6 around the phosphorylation site of Y219.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A