

CCNE2 Antibody



PACO05695

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CCNE2

Uniprot

O96020

Synonyms:

CCNE2; G1/S-specific cyclin-E2

Immunogen:

Synthesized peptide derived from human Cyclin E2 around the non-phosphorylation site of T392.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A