

DYRK1B Antibody



PACO05700

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

DYRK1B

Uniprot

Q9Y463

Synonyms:

DYRK1B; MIRK; Dual specificity tyrosine-phosphorylation-regulated kinase 1B;
Minibrain-related kinase; Mirk protein kinase

Immunogen:

Synthesized peptide derived from the Internal region of human Dyrk1B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A