

Phospho-RPS6KA5 (S360) Antibody



PACO05728

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

RPS6KA5

Uniprot

O75582

Synonyms:

RPS6KA5; MSK1; Ribosomal protein S6 kinase alpha-5; S6K-alpha-5; 90 kDa ribosomal protein S6 kinase 5; Nuclear mitogen- and stress-activated protein kinase 1; RSK-like protein kinase; RSKL

Immunogen:

Synthesized peptide derived from human MSK1 around the phosphorylation site of S360.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A