

CCDC102A Antibody



PACO05737

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

CCDC102A

Uniprot

Q96A19

Synonyms:

CCDC102A; Coiled-coil domain-containing protein 102A

Immunogen:

Synthesized peptide derived from the Internal region of human CCDC102A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A