

RPL23 Antibody



PACO05748

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

RPL23

Uniprot

P62829

Synonyms:

RPL23; 60S ribosomal protein L23; 60S ribosomal protein L17

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L23.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A