

# SLC22A6 Antibody



PACO05749

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

SLC22A6

**Uniprot**

O95742

**Synonyms:**

Solute carrier family 22 member 6; Organic anion transporter 1; Hoat1; Renal organic anion transporter 1; hROAT1; PAH transporter; hPAHT; SLC22A6; OAT1; PAHT

**Immunogen:**

Synthesized peptide derived from the Internal region of human OAT1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A