

# RAB33A Antibody



PACO05750

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

RAB33A

**Uniprot**

Q14088

**Synonyms:**

RAB33A; RABS10; Ras-related protein Rab-33A; Small GTP-binding protein S10

**Immunogen:**

Synthesized peptide derived from the Internal region of human Rab 33A.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A