

COX5A Antibody



PACO05752

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

COX5A

Uniprot

P20674

Synonyms:

COX5A; Cytochrome c oxidase subunit 5A; mitochondrial; Cytochrome c oxidase polypeptide Va

Immunogen:

Synthesized peptide derived from the Internal region of human COX5a.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A