

SYTL4 Antibody



PACO05757

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SYTL4

Uniprot

Q96C24

Synonyms:

SYTL4; Synaptotagmin-like protein 4; Exophilin-2; Granuphilin

Immunogen:

Synthesized peptide derived from the Internal region of human Granuphilin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A