

MP68 Antibody



PACO05778

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

MP68

Uniprot

P56378

Synonyms:

MP68; C14orf2; 6.8 kDa mitochondrial proteolipid

Immunogen:

Synthesized peptide derived from the Internal region of human MP68.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A