

RIN1 Antibody



PACO05790

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

RIN1

Uniprot

Q13671

Synonyms:

RIN1; Ras and Rab interactor 1; Ras inhibitor JC99; Ras interaction/interference protein 1

Immunogen:

Synthesized peptide derived from the C-terminal region of human RIN1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A