

NCOR1 Antibody



PACO05791

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

NCOR1

Uniprot

O75376

Synonyms:

NCOR1; KIAA1047; Nuclear receptor corepressor 1; N-CoR; N-CoR1

Immunogen:

Synthesized peptide derived from the N-terminal region of human N-CoR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A