

PACO05794

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:10000, IF:1:200-1:1000

Protein Background:**Gene ID:**

CELSR1

Uniprot

Q9NYQ6

Synonyms:

CELSR1; CDHF9; FMI2; Cadherin EGF LAG seven-pass G-type receptor 1; Cadherin family member 9; Flamingo homolog 2; hFmi2

Immunogen:

Synthesized peptide derived from the Internal region of human CDHF9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A