

GNA14 Antibody



PACO05795

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

GNA14

Uniprot

O95837

Synonyms:

GNA14; Guanine nucleotide-binding protein subunit alpha-14; G alpha-14; G-protein subunit alpha-14

Immunogen:

Synthesized peptide derived from the N-terminal region of human G α 14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A