

# USP6NL Antibody



PACO05805

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

USP6NL

**Uniprot**

Q92738

**Synonyms:**

USP6NL; KIAA0019; USP6 N-terminal-like protein; Related to the N-terminus of tre; RN-tre

**Immunogen:**

Synthesized peptide derived from the Internal region of human RN-tre.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A