

ACYP2 Antibody



PACO05816

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ACYP2

Uniprot

P14621

Synonyms:

ACYP2; ACYP; Acylphosphatase-2; Acylphosphatase; muscle type isozyme; Acylphosphate phosphohydrolase 2

Immunogen:

Synthesized peptide derived from the Internal region of human ACYP2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A