

CEP350 Antibody



PACO05829

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

CEP350

Uniprot

Q5VT06

Synonyms:

CEP350; CAP350; KIAA0480; GM133; Centrosome-associated protein 350; Cep350; Centrosome-associated protein of 350 kDa

Immunogen:

Synthesized peptide derived from the Internal region of human CEP350.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A