

# REC8 Antibody



PACO05862

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

REC8

**Uniprot**

O95072

**Synonyms:**

REC8; REC8L1; Meiotic recombination protein REC8 homolog; Cohesin Rec8p

**Immunogen:**

Synthesized peptide derived from the Internal region of human Rec8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A