

STAT2 Antibody



PACO05866

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

WB:1:1000-1:2000

Protein Background:

Gene ID:

STAT2

Uniprot

P52630

Synonyms:

Signal transducer and activator of transcription 2

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A