

CDH4 Antibody



PACO05892

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

CDH4

Uniprot

P55283

Synonyms:

CDH4; Cadherin-4; Retinal cadherin; R-CAD; R-cadherin

Immunogen:

Synthesized peptide derived from the C-terminal region of human R-cadherin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A