

IRS1 Antibody



PACO05917

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

IRS1

Uniprot

P35568

Synonyms:

IRS1; Insulin receptor substrate 1; IRS-1

Immunogen:

Synthesized peptide derived from human IRS-1 around the non-phosphorylation site of S323.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A