

ITGB2 Antibody



PACO05933

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:5000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ITGB2

Uniprot

P05107

Synonyms:

ITGB2; CD18; MFI7; Integrin beta-2; Cell surface adhesion glycoproteins LFA-1/CR3/p150; 95 subunit beta; Complement receptor C3 subunit beta; CD antigen CD18

Immunogen:

Synthesized peptide derived from human Integrin β 2 around the non-phosphorylation site of T758.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A