

JAK2 Antibody



PACO05948

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

JAK2

Uniprot

O60674

Synonyms:

JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2

Immunogen:

Synthesized peptide derived from human JAK2 around the non-phosphorylation site of Y221.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A