

# Phospho-JUNB (S79) Antibody



PACO05953

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IP

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300, IP:2-5ug

**Protein Background:****Gene ID:**

JUNB

**Uniprot**

P17275

**Synonyms:**

JUNB; Transcription factor jun-B

**Immunogen:**

Synthesized peptide derived from human Jun B around the phosphorylation site of S79.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A