

Phospho-KCNB1 (S567) Antibody



PACO05960

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

KCNB1

Uniprot

Q14721

Synonyms:

KCNB1; Potassium voltage-gated channel subfamily B member 1; Delayed rectifier potassium channel 1; DRK1; h-DRK1; Voltage-gated potassium channel subunit Kv2.1

Immunogen:

Synthesized peptide derived from human Kv2.1 around the phosphorylation site of S567.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A