

# SFT2D2 Antibody



PACO05964

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

SFT2D2

**Uniprot**

O95562

**Synonyms:**

SFT2D2; Vesicle transport protein SFT2B; SFT2 domain-containing protein 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human SFT2D2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A