

# OR56B1 Antibody



PACO05976

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:10000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

OR56B1

**Uniprot**

Q8NGI3

**Synonyms:**

OR56B1; OR56B1P; Olfactory receptor 56B1; Olfactory receptor OR11-65

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 56B1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A