

PAK2 Antibody



PACO05981

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

PAK2

Uniprot

Q13177

Synonyms:

PAK2; Serine/threonine-protein kinase PAK 2; gamma-PAK; PAK65; S6/H4 kinase; p21-activated kinase 2; PAK-2; p58

Immunogen:

Synthesized peptide derived from human PAK γ around the non-phosphorylation site of S141.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A