

CXCL10 Antibody



PACO06008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

CXCL10

Uniprot

P02778

Synonyms:

CXCL10; INP10; SCYB10; C-X-C motif chemokine 10; 10 kDa interferon gamma-induced protein; gamma-IP10; IP-10; Small-inducible cytokine B10

Immunogen:

Synthesized peptide derived from the Internal region of human IP-10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A