

OR51B6 Antibody



PACO06011

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR51B6

Uniprot

Q9H340

Synonyms:

OR51B6; Olfactory receptor 51B6; Odorant receptor HOR5'beta6

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 51B6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A