

# TGIF1 Antibody



PACO06023

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

TGIF1

**Uniprot**

Q15583

**Synonyms:**

TGIF1; TGIF; Homeobox protein TGIF1; 5'-TG-3'-interacting factor 1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human TGIF.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A